CEKA SITE

3 °C
37 °F

9 °C
48 °F

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EN

Two-component composite for adhering prosthetic attachments.

Directions for use

CEKA SITE is used to adhere metal components directly on the working model, ensuring an absolute passive fit. A high proportion of inorganic fillers guarantees a high final hardness producing the stability required for laboratory preparation.

Indications

- Fixing telescopic or conical crowns to chrome-cobalt frameworks
- Fixing prefabricated attachments to chrome-cobalt frameworks

Constructing the joint

The areas to be adhered must be designed so as to be retentive and strong. When the secondary component is placed on the primary component, an even, thin gap must be present. This gap, which ensures that a passive fit is obtained, is then filled with CEKA SITE. When adhering prefabricated attachments, undercuts must be created in the joint areas and the manufacturer's instructions observed.

Preparing the bonding surfaces

After ensuring that both the primary and secondary components fit passively, create undercuts in them and sandblast with an aluminium oxide (max. 125 µm). After steam cleaning, condition the metal surfaces with a metal bonding system.

Dosage

After putting on the mixing tip and pushing the plunger of the syringe, the two components get automatically mixed. Please leave the mixing tip on the syringe as a closure after use.

Working time after starting mixing: approx. 3 min.

Setting commences after approx. 4½ min.

Setting completed after approx. 8 min.

These times assume a room temperature of 22 °C.

Application

Use a spatula to apply CEKA SITE to the female and male. The metal components should be repositioned on the model and fixed into place. Ensure that the areas to be bonded are completely filled with CEKA SITE and that there are not any air bubbles or contamination.

Important

The sections being joined must not be moved until CEKA SITE has set completely.

Please note

- CEKA SITE is dimensionally stable up to 120 °C.
- CEKA SITE is easy to grind.

Storage

If refrigerated at 3-9 °C, this product can be kept for 2 years. Do not use after the expiry date.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects - even if it is doubtful by our product - please contact us.

Contra-indications / Interactions

If a patient has known allergies against or hypersensitiveness towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross reactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above information if the medical device is used to produce a special model. Please also consider the safety data sheets. Avoid contact of unpolymerized material with skin and mucous membranes.

Troubleshooting

Working time is too short.	High ambient temperatures caused by e.g. summer heat or central heating will decrease the working time.
Working time is too long.	If you use the product directly from the refrigerator, this will increase the working time.
Material does not polymerize sufficiently.	- Substances containing eugenol/oil of wintergreen inhibit polymerization. - The material was not intermixed properly. Only use original mixing tips.
Material has already polymerized.	After the last use, the cap was put back onto the syringe instead of leaving the mixing tip on it (contamination of components A and B).
Construction does not fit.	- The objects were moved before the polymerization was completed. - Material was already too hard: Please bond fewer parts at a time.

Contents / Re-ordering

CEKA SITE B (9795832) 2 x 2 g composite, 6 automix tips

6705 (9795437) CEKA Multi automix tips 20 pcs

Hazard

Contains tetramethylene dimethacrylate, dibenzoyl peroxide.

Warning

May cause an allergic skin reaction. Avoid breathing vapours/spray. Wear protective gloves. If skin irritation or rash occurs: Get medical advice attention.

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