

FABRICATION OF A NEW PROSTHESIS WITH PRECI-BAR

Article numbers:	1100/H/MR/OR	1100/H/MR/PA	1100/H/MR/IN
	1102/H/MR/OR	1102/H/MR/PA	1102/H/MR/IN
	1103/B/OR	1103/B/PA	1012/B
	1103/BL/OR	1103/BL/PA	
	1101/B/OR	1101/B/PA	1013/B
	1101/BL/OR	1101/BL/PA	
	1100/B/OR	1100/B/PA	1011/B
	1100/BL/OR	1100/BL/PA	

1. Take a full arch impression.
2. Pour a stone model. The entire bar area must be made of epoxy material or self-curing acrylic resin.
3. Reduce the corresponding rider to the required length (according to the bar) **(1100/H/MR/OR, 1100/H/MR/PA, 1100/H/MR/IN, 1102/H/MR/OR, 1102/H/MR/PA, 1102/H/MR/IN)**.
4. Place the rider on the model. When processing a resilient PRECI-BAR, the space maintainer must be used.
5. Prepare the wax setup for investing.
6. Separate the reproduced bar from the acrylic resin. Complete the prosthesis and remove the space maintainer if used.
7. Clean the rider and activate it if necessary (see INFO 108).

Side effects, warnings and precautions

- The attachments are intended for single use.
- The products are non-sterile.
- There is a risk of poor fit when patient conditions change.
- Bacterial adhesion can be avoided by applying hygiene measures.
- Inappropriate use or bad manufacturing can lead to premature wear of the attachments.
- The functionality of the attachments will be adversely affected by traumas such as grinding and bruxism.
- For the purpose of traceability we advise you to record the lot number of the applied products in the patient file.
- Do not heat items containing titanium.
- Do not use items containing nickel in case of nickel allergy.
- The accessories RE H 79 and H 35 must be used outside the mouth.