

PROCESSING THE MALES OF PRECI-CLIX TI-P/P-P

1. Determine the path of insertion of the prosthesis.
2. Select the appropriate **RE 04 TI** plastic male keeper and press the **1281 C** male into the cavity.
3. Place the **1201 P** paralleling mandrel carefully on top of the male.
4. Adapt the plastic pattern according to the shape, position and height of the abutment crowns.
5. Fix the male keeper with wax to the abutment crowns. Remove the paralleling mandrel.
6. Sprue the crowns and prepare for investing.
7. Invest, burnout, preheat carefully, and cast in a hard alloy. Do not use fast heating investments.
8. Sandblast the casting without damaging the male. Use the **1230** cup burr.
9. Carefully finish the casting. Do not remove too much material from the ball. Thoroughly polish.

Side effects, warnings and precautions

- The attachments are intended for single use.
- The products are non-sterile.
- There is a risk of poor fit when patient conditions change.
- Bacterial adhesion can be avoided by applying hygiene measures.
- Inappropriate use or bad manufacturing can lead to premature wear of the attachments.
- The functionality of the attachments will be adversely affected by traumas such as grinding and bruxism.
- For the purpose of traceability we advise you to record the lot number of the applied products in the patient file.
- Do not heat items containing titanium.
- Do not use items containing nickel in case of nickel allergy.
- The accessories RE H 79 and H 35 must be used outside the mouth.