3C-BOND

3C-BOND
140-3CB-20

20 g

CEKA AXIAL
PRECI-CLIX AXIAL
PRECI-CLIX RADICULAR
PRECI-BALL

CEKA EXTRACORONAL
PRECI-VERTIX
PRECI-CLIX EC
PRECI-SAGIX
PRECI-52

PRECI-BAR
PRECI-HORIX
PRECI-CLIP

PRECI-PROFILE
PRECI-POST
CEKA SOL
CEKA SITE
CEKA BOND
3C-BOND
PERMA-RET
PRECI-SEP
PLASTICWAX
EXPANDO
CEKA Multi (O)
MEASURING GAUGE
CKPL DIGITAL LIBRARY

01.
AXIAL
ATTACHMENTS

02.
EXTRACORONAL
ATTACHMENTS

03.
PROFILE
ATTACHMENTS

SPECIALTIES

SPECIALTIES

3C-BOND EN-DE-FR

WWW.CKPL.EU

SIDE EFFECTS, WARNINGS AND PRECAUTIONS

The attachments are intended for single use.
The products are non-sterile.
There is a risk of poor fit when patient conditions change.
Bacterial adhesion can be avoided by applying hygiene measures.
Inappropriate use or bad manufacturing can lead to premature wear of the attachments.
The functionality of the attachments will be adversely affected by traumas such as grinding and bruxism.
For the purpose of traceability we advise you to record the lot number of the applied products in the patient file.
Do not heat products containing titanium.
Do not use products containing nickel in case of nickel allergy.
The accessories RE H 79 and H 35 must be used outside the mouth.

CEKA
ATTACHMENTS
PRECI-LINE

Prospectus no 275.M1-ed.F
© Copyright ALPHADENT NV 09/2017 - BTW BE 0414.789.321 RPR Kortrijk
ALPHADENT NV • Mannebeekstraat 33, 8790 Waregem, Belgium
T +32 56 629 531 - E info@ckpl.eu

CEKA
ATTACHMENTS
PRECI-LINE

3C-BOND

- ✓ For a warm and natural porcelain shade
- ✓ Exceptional bonding
- ✓ No greening

3C-BOND is a bonding material for porcelain to - titanium-free - non-precious and chrome-cobalt alloys. 3C-BOND avoids the formation of green oxides and improves the bonding with the metal. Use with conventional porcelain.

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INSTRUCTIONS FOR USE

Sandblast with aluminium oxide.
Clean in an ultrasonic bath or by boiling for 10 minutes in distilled water, or steam clean.
The metal structure must be completely dry and should only be touched with clean tweezers.
Mix 3C-BOND with distilled water and apply it in a thin layer. Cover also marginal areas.

Never re-use already mixed 3C-BOND!

Fire 3C-BOND under vacuum in the same way as opaque.
Preheating temperature: 650 °C (1200 °F)
Temperature rise: 55 °C (130 °F) per minute
Final temperature: 980 °C (1800 °F)
After firing, immediately remove the metal structure from the furnace and bench cool.
3C-BOND should have a yellowish, slightly shiny surface.
Continue with the opaque firing cycle as recommended by the porcelain manufacturer.

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