


CEKA SOL

FL INCORPORATED
EINGEARBEITET FL MI EL
INGE ERK LOEIMIDDEL



CEKA AXIAL
PRECI-CLIP AXIAL
PRECI-CLIP RADICULAR
PRECI-BALL

CEKA EXTRACORONAL
PRECI-CLIP
PRECI-CLIP EC
PRECI-CLIP AGI
PRECI-52

PRECI-BAR
PRECI-HORICON
PRECI-CLIP

PRECI-PROFILE
PRECI-PO

CEKA SOL
CEKA TIE
CEKA BOND
3C-BOND
PERMA-RE
PRECI-EP
PLAIC A
E PANDO
CEKA M (O)
MEARING GAUGE
CKPL DIGITAL LIBRARY

01.
AXIAL
ATTACHMENTS

02.
EXTRACORONAL
ATTACHMENTS

03.
PROFILE
ATTACHMENTS

SPECIALTIES

SPECIALTIES

CEKA SOL EN-DE-NL

WWW.CKPL.EU

SIDE EFFECTS, WARNINGS AND PRECAUTIONS

The attachments are intended for single use.
The products are non-sterile.
There is a risk of poor fit when patient conditions change.
Bacterial adhesion can be avoided by applying hygiene measures.
Inappropriate use or bad manufacturing can lead to premature wear of the attachments.
The functionality of the attachments will be adversely affected by traumas such as grinding and bruxism.
For the purpose of traceability we advise you to record the lot number of the applied products in the patient file.
Do not heat products containing titanium.
Do not use products containing nickel in case of nickel allergy.
The accessories RE H 79 and H 35 must be used outside the mouth.



CEKA
ATTACHMENTS
PRECI-LINE

SWISS
MADE

CE0197



Prospectus no 268.M1-ed.D
C ALPHADEN N 03/2017 - B BE 0414.789.321 RPR K
ALPHADEN N M 33, 8790 ,B
+32 56 629 531 - E @ .

CEKA
ATTACHMENTS
PRECI-LINE

LOW-FUSING TUBE SOLDER WITH FLUX INCORPORATED



A applications. The incorporated flux provides for



Simple & easy!

		PRECIOUS	NON-PRECIOUS	CR-CO
CEKA SOL	# CS00	X	X	X
780-820 C (1436-1508 F) A 80 - C 13 - I 5 - 2				
CEKA SOL FILIGRAN	# CS00 F	X	X	X
780-820 C (1436-1508 °F) (contains less flux) A 80 - C 13 - I 5 - 2				
CEKA SOL W	# CS22	X	X	X
865-930 C (1589-1706 F) A 76,5 - C 0,4 - 6,5 - N 16,6				
CEKA SOL L	# CS21	X		X
615-630 C (1139-1166 F) A 42 - A 38,5 - 11 - I 8,5				
CEKA SOL G	# CS20	X		
755-820 C (1391-1508 F) A 68 - A 13,3 - C 12,8 - 1 - 4,9				
CEKA SOL CR	# CS23			X
1000-1030 C (1832-1886 F) A 10 - A 37 - P 37 - C 12 - 4				

ADDITION ALLOY FOR LASER WELDING WITHOUT FLUX

CEKA SOL T1 # CS31
Round 8 x 25 cm Diameter: 0.4 mm

NIEDRIGSCHMELZENDES LOT IN RÖHRENFORM, MIT EINGEARBEITETEM FLUSSMITTEL



Ein komplettes Angebot von niedrigschmelzenden Loten für verschiedene Anwendungen. Durch das eingearbeitete Fluxmittel werden die Lötstellen fest miteinander verbunden.



Einfach problemlos!

		EDELMETALL	NICHTEDELMETALL	CR-CO
CEKA SOL	# CS00	X	X	X
G 780-820 C A 80 - C 13 - I 5 - 2				
CEKA SOL FILIGRAN	# CS00 F	X	X	X
G 780-820 C () A 80 - C 13 - I 5 - 2				
CEKA SOL W	# CS22	X	X	X
865-930 C A 76,5 - C 0,4 - 6,5 - N 16,6				
CEKA SOL L	# CS21	X		X
G 615-630 C A 42 - A 38,5 - 11 - I 8,5				
CEKA SOL G	# CS20	X		
G 755-820 C A 68 - A 13,3 - C 12,8 - 1 - 4,9				
CEKA SOL CR	# CS23			X
1000-1030 C A 10 - A 37 - P 37 - C 12 - 4				

ADDITIONSLEGIERUNG FÜR LASERSCHWEISSEN OHNE FLUSSMITTEL

CEKA SOL T1 # CS31
Rund 8 x 25 cm Durchmesser: 0,4 mm

LAAGSMELTEND SOLDEER IN STAAFVORM, MET INGEWERKT VLOEIMIDDEL



Een compleet palet van laagsmeltende solderen voor diverse toepassingen. Het ingewerkte vloeimiddel garandeert een vlot solderen van verschillende types



Eenvoudigweg probleemloos!

		EDELMETAAL	NIET-EDELMETAAL	CR-CO
CEKA SOL	# CS00	X	X	X
G 780-820 C A 80 - C 13 - I 5 - 2				
CEKA SOL FILIGRAN	# CS00 F	X	X	X
G 780-820 C (bevat minder vloeimiddel) A 80 - C 13 - I 5 - 2				
CEKA SOL W	# CS22	X	X	X
865-930 C A 76,5 - C 0,4 - 6,5 - N 16,6				
CEKA SOL L	# CS21	X		X
G 615-630 C A 42 - A 38,5 - 11 - I 8,5				
CEKA SOL G	# CS20	X		
G 755-820 C A 68 - A 13,3 - C 12,8 - 1 - 4,9				
CEKA SOL CR	# CS23			X
1000-1030 C A 10 - A 37 - P 37 - C 12 - 4				

ADDITIONELEGERING VOOR LASERLASSEN ZONDER VLOEIMIDDEL

CEKA SOL T1 # CS31
Wit, zuiver titaan
Rond 8 x 25 cm Diameter: 0,4 mm