

PRECI-POST M



PRECI-POST CC



CEKA AXIAL
PRECI-CLIX AXIAL
PRECI-CLIX RADICULAR
PRECI-BALL

CEKA EXTRACORONAL
PRECI-VERTIX
PRECI-CLIX EC
PRECI-SAGIX
PRECI-52

PRECI-BAR
PRECI-HORIX
PRECI-CLIP

PRECI-PROFILE
PRECI-POST
CEKA SOL
CEKA SITE
CEKA BOND
3C-BOND
PERMA-RET
PRECI-SEP
PLASTICWAX
EXPANDO
CEKA Multi (O)
MEASURING GAUGE
CKPL DIGITAL LIBRARY

01.
AXIAL
ATTACHMENTS

02.
EXTRACORONAL
ATTACHMENTS

03.
PROFILE
ATTACHMENTS

SPECIALTIES

SPECIALTIES

PRECI-POST ENGLISH

WWW.CKPL.EU

SIDE EFFECTS, WARNINGS AND PRECAUTIONS

The attachments are intended for single use.
The products are non-sterile.
There is a risk of poor fit when patient conditions change.
Bacterial adhesion can be avoided by applying hygiene measures.
Inappropriate use or bad manufacturing can lead to premature wear of the attachments.
The functionality of the attachments will be adversely affected by traumas such as grinding and bruxism.
For the purpose of traceability we advise you to record the lot number of the applied products in the patient file.
Do not heat products containing titanium.
Do not use products containing nickel in case of nickel allergy.
The accessories RE H 79 and H 35 must be used outside the mouth.



CEKA
ATTACHMENTS
PRECI-LINE

SWISS
MADE

CE0197



Prospectus no 261.EN-ed.H
© Copyright ALPHADENT NV 03/2017 - BTW BE 0414.789.321 RPR Kortrijk
ALPHADENT NV • Mannebeekstraat 33, 8790 Waregem, Belgium
T +32 56 629 531 - E info@ckpl.eu

CEKA
ATTACHMENTS
PRECI-LINE

PRECI-POST

PRECI-POST M

Rational system with a conicity of 2°30' over the entire useful length. Form and dimensions according to Dr Mooser. For post-copings and temporary post cores with one or more roots. Yellow and blue colour codification.



PRECI-POST M SET



PRECI-POST M LABORATORY POST



PRECI-POST M IMPRESSION POST

PRECI-POST M



PRECI-POST M SET # 2001

- 1 yellow reamer
- 1 blue reamer
- 25 yellow opaque impression posts
- 25 blue opaque impression posts
- 25 yellow transparent laboratory posts
- 25 blue transparent laboratory posts

PRECI-POST M LABORATORY POST

Transparent plastic posts for lab use only.

- # 2003GL 50 pcs yellow tip \varnothing 0.80 mm
cone length 10.70 mm, total length 18.60 mm
- # 2003BL 50 pcs blue tip \varnothing 1.00 mm
cone length 12.30 mm, total length 22.90 mm
- # 2003L 2 x 50 pcs combination package
yellow & blue

PRECI-POST M IMPRESSION POST

Opaque plastic impression posts (with retention), to be used with any impression technique (including silverplating and epoxy)

- # 2002G 50 pcs yellow tip \varnothing 0.80 mm
length 14.10 mm
- # 2002B 50 pcs blue tip \varnothing 1.00 mm
length 16.70 mm
- # 2002 2 x 50 pcs combination package
yellow & blue

PRECI-POST M REAMER

Reamer according to the Mooser system

- # 2004G 1 pc yellow tip \varnothing 0.80 mm
- # 2004B 1 pc blue tip \varnothing 1.00 mm

PRECI-POST CC



PRECI-POST CC SET

For lab use only.

- # 2050 50 yellow posts, 50 red posts,
50 blue posts, 50 green posts
- # 2060 20 yellow posts, 15 red posts,
10 blue posts, 5 green posts

PRECI-POST CC PLASTIC POST

Plastic posts for lab use only.

- # 2051A 50 pcs yellow \varnothing 1.40 mm, tip \varnothing 0.65 mm
cone length 4.60 mm, total length 16.60 mm
- # 2052B 50 pcs red \varnothing 1.60 mm, tip \varnothing 0.70 mm
cone length 5.60 mm, total length 18.60 mm
- # 2053C 50 pcs blue \varnothing 1.80 mm, tip \varnothing 0.75 mm
cone length 6.40 mm, total length 20.00 mm
- # 2054D 50 pcs green \varnothing 2.00 mm, tip \varnothing 0.80 mm
cone length 7.05 mm, total length 22.40 mm

PRECI-POST CC REAMER

Cylindro-conical reamer

- # 2951A 1 pc yellow tip \varnothing 0.65 mm
- # 2952B 1 pc red tip \varnothing 0.70 mm
- # 2953C 1 pc blue tip \varnothing 0.75 mm
- # 2954D 1 pc green tip \varnothing 0.80 mm

PRECI-POST CC

Cylindro-conical system: 5°30' at the base and then cylindrical over the entire useful length for optimum retention. For post-copings and temporary post cores with one or more roots. Yellow, red, blue, and green colour codification.



PRECI-POST CC SET



PRECI-POST CC REAMER